

# NOVA

Advanced Ultra Short Progressive

# Petite



Sometimes fashion frames are too small for regular progressive lenses and the wearer is forced to select a bigger frame. Nova-Petite, specially developed for small frames, offers more options and flexibility to select the ideal progression length perfectly adapted to each frame.

**NOVA** OPTICAL LABS  
TRULY ALBERTAN INDEPENDENT

Calgary 1-403-272-2007

Red Deer 1-403-346-0999

Lethbridge 1-403-329-0041

[www.Novaopticalabs.ca](http://www.Novaopticalabs.ca)

# NOVA-PETITE

## Advanced Ultra Short

Ideal as a daily all-purpose progressive lens and for those wearers who want to wear very small frames



Ideal solution for users with demanding requirements for large visual fields at any distance. These patients are characterized by asking for a comfortable lens while at the same time wanting excellent vision of objects located at any distance.

### Advantages

Developed for patients who want to wear small frames

Fast transition between far and near vision

Small progression length to adaptable to the smallest frames

Available in four progression lengths

High precision and high personalization due to Digital Ray-Path Technology

Clear vision in every gaze direction

Oblique astigmatism minimized


Variable inset: Automatic and manual

Frame shape optimization available.

Vertex distance	✓
Near working distance	✓
Pantoscopic angle	✓
Wrapping angle	✓
IPD	✓
SEGHT	✓
HBOX	✓
VBOX	✓

Minimum Fitting Height MFH
10mm
11mm
12mm
13mm

### Design Performance



Performance at each distance	
Far	★★★★☆☆
Intermediate	★★★★☆☆
Near	★★★★☆☆

**NOVA** OPTICAL LABS  
TRULY ALBERTAN INDEPENDENT